

GOAT ANTI-URP2 / KINDLIN-3 ANTIBODY

SKU: EB09539



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	FERMT3 fermitin family homolog 3 (Drosophila) KIND3 MGC10966 MIG-2 MIG2B UNC112C URP2 URP2SF UNC-112 related protein 2 fermitin family homolog 3 kindlin 3 kindlin-3
Accession ID	NP_113659.3; NP_848537.1
Blocking Peptide	EBP09539
Immunogen	Peptide with sequence C-ERARGEELDED, from the C Terminus of the protein sequence according to NP_113659.3; NP_848537.1.
Product Comments	This antibody is expected to recognize both isoforms (NP_113659.3, NP_848537.1).
Peptide Sequence	C-ERARGEELDED
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog
Reactive Species	Human
Human Gene ID	83706
Mouse Gene ID	108101
Rat Gene ID	309186
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:2000.
Western Blot	Approx 75kDa band observed in Human Peripheral Blood Mononucleocytes lysates (calculated MW of 75.4kDa according to NP_113659.3). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

