

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-URP2 / KINDLIN-3 ANTIBODY

SKU: EB09539



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / FERMT3|fermitin family homolog 3

(Drosophila)|KIND3|MGC10966|MIG-2|MIG2B|UNC112C|URP2|URP2SF|UNC-112 related protein 2|fermitin family **Alias**

homolog 3|kindlin 3|kindlin-3 Names

Accession

NP_113659.3; NP_848537.1

Blocking

ID

EBP09539 **Peptide**

Peptide with sequence C-ERARGEELDED, from the C Terminus of the protein sequence according to **Immunogen**

NP_113659.3; NP_848537.1.

Product Comments

This antibody is expected to recognize both isoforms (NP 113659.3, NP 848537.1).

Peptide

Sequence

C-ERARGEELDED

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Mouse, Rat, Dog

Species Reactive

Human **Species**

Human

83706

Gene ID

Mouse 108101 Gene ID **Rat Gene ID** 309186

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Grade **ELISA**

Detection Antibody detection limit dilution 1:2000.

Limit

Western Approx 75kDa band observed in Human Peripheral Blood Mononucleocytes lysates (calculated MW of 75.4kDa

according to NP_113659.3). Recommended concentration: 1-3µg/ml. **Blot**

Application

Type

Pep-ELISA, WB



 $\textbf{Email:} \ \underline{\textbf{customerservice@vectorlabs.com}}$

Telephone: (650) 697-3600

GALLERY IMAGES

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa