

GOAT ANTI-USP11 / UHX1 ANTIBODY

SKU: EB05984



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	USP11 UHX1 ubiquitin specific protease 11 deubiquitinating enzyme 11 ubiquitin thiolesterase 11 ubiquitin carboxyl-terminal hydrolase, X-linked ubiquitin specific peptidase 11 RP4-659F15.2 ubiquitin-specific processing protease 11
Accession ID	NP_004642.2
Blocking Peptide	EBP05984
Immunogen	Peptide with sequence CSSPPSSEFMDVN, from the C Terminus of the protein sequence according to NP_004642.2.
Peptide Sequence	CSSPPSSEFMDVN
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	8237
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Western Blot	Approx 120kDa band observed in lysates of cell line Jurkat (calculated MW of 110kDa according to NP_004642.2). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Ideguchi et al, Biochem J. 2002 Oct 1;367(Pt 1):87-95.; PMID: 12084015). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

