

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-USP11 / UHX1 ANTIBODY

SKU: EB05984



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa



Telephone: (650) 697-3600



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

USP11|UHX1|ubiquitin specific protease 11|deubiquitinating enzyme 11|ubiquitin thiolesterase 11|ubiquitin carboxyl-terminal hydrolase, X-linked|ubiquitin specific peptidase 11|RP4-659F15.2|ubiquitin-specific

Alias Names processing protease 11

Accession

ID

NP_004642.2

Blocking

Peptide

EBP05984

Immunogen

Peptide with sequence CSSPPSSEFMDVN, from the C Terminus of the protein sequence according to

NP 004642.2.

Peptide

Sequence

CSSPPSSEFMDVN

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Instructions

Refrigerated

Predicted

Species

Human

Human

8237

Reactive

Species

Human

Gene ID

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA

Detection

Western

Antibody detection limit dilution 1:8000.

Limit

Approx 120kDa band observed in lysates of cell line Jurkat (calculated MW of 110kDa according to NP_004642.2). The observed molecular weight corresponds to earlier findings in literature with different

Blot antibodies (Ideguchi et al, Biochem J. 2002 Oct 1;367(Pt 1):87-95.; PMID: 12084015). Recommended

concentration: 1-3µg/ml.

Application

Type

Pep-ELISA, WB

GALLERY IMAGES





Telephone: (650) 697-3600



250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa