

# GOAT ANTI-USP16 ANTIBODY

**SKU:** EB08657



## SPECIFICATIONS

<b>Formulation</b>	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
<b>Unit Size</b>	100 µg
<b>Storage Instructions</b>	Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Synonym / Alias Names</b>	ubiquitin-specific processing protease 16 ubiquitin thiolesterase 16 ubiquitin specific protease 16 ubiquitin carboxyl-terminal hydrolase 16 deubiquitinating enzyme 16 UBP-M ubiquitin specific peptidase 16 USP16
<b>Accession ID</b>	NP_001027582.1; NP_001001992.1
<b>Blocking Peptide</b>	EBP08657
<b>Immunogen</b>	Peptide with sequence NVEWNICQDCKTDNK, from the internal region of the protein sequence according to NP_001027582.1; NP_001001992.1.
<b>Product Comments</b>	Reported variants represent identical protein: NP_001027582.1 NP_006438.1
<b>Peptide Sequence</b>	NVEWNICQDCKTDNK
<b>Purification Method</b>	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
<b>Shipping Instructions</b>	Refrigerated
<b>Predicted Species</b>	Human, Mouse, Rat, Horse, Dog, Cow
<b>Human Gene ID</b>	10600
<b>Mouse Gene ID</b>	74112
<b>Rat Gene ID</b>	288306
<b>Product Grade</b>	<a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png</a>
<b>ELISA Detection Limit</b>	Antibody detection limit dilution 1:32000.
<b>Application Type</b>	Pep-ELISA

## DOCUMENTS

- [Data Sheet](#)