Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-USP16 ANTIBODY

SKU: EB08657

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias

ubiquitin-specific processing protease 16|ubiquitin thiolesterase 16|ubiquitin specific protease 16|ubiquitin carboxyl-terminal hydrolase 16|deubiquitinating enzyme 16|UBP-M|ubiquitin specific peptidase 16|USP16

Names Accession

NP_001027582.1; NP_001001992.1

Blocking

Peptide

Peptide with sequence NVEWNICQDCKTDNK, from the internal region of the protein sequence according to Immunogen

NP 001027582.1; NP 001001992.1.

Product

Reported variants represent identical protein: NP_001027582.1 NP_006438.1 Comments

Peptide

Sequence

NVEWNICODCKTDNK

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

EBP08657

Shipping

Instructions

Refrigerated

Predicted

Species

Human, Mouse, Rat, Horse, Dog, Cow

Human **Gene ID**

10600

Mouse

74112

Gene ID

Rat Gene ID 288306

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

Grade **ELISA**

Detection

Antibody detection limit dilution 1:32000.

Limit

Application Pep-ELISA

Type





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600