

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-USP18 / UBP43 ANTIBODY

SKU: EB05996

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias UBP43|ISG43|ubiquitin specific protease 18|USP18

Names

Accession NP 059110

Immunogen Peptide with sequence SKAFGLLRQICQS, from the N Terminus of the protein sequence according to NP_059110.

Peptide SKAFGLLRQICQS

Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Human 11274

Gene ID

Product https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:2000.

Limit

Blot

Western Blot: No signal obtained yet but low background observed in Human Liver extracts at up to 2µg/ml. We Western

would appreciate any feedback from people in the field - have any results been reported with other

antibodies/lysates?