



GOAT ANTI-USP18 / UBP43 ANTIBODY

SKU: EB05996



SPECIFICATIONS

| | |
|------------------------------|---|
| Formulation | Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. |
| Unit Size | 100 µg |
| Storage Instructions | Aliquot and store at -20°C. Minimize freezing and thawing. |
| Synonym / Alias Names | UBP43 ISG43 ubiquitin specific protease 18 USP18 |
| Accession ID | NP_059110 |
| Immunogen | Peptide with sequence SKAFGLLRQICQS, from the N Terminus of the protein sequence according to NP_059110. |
| Peptide Sequence | SKAFGLLRQICQS |
| Purification Method | Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide. |
| Shipping Instructions | Refrigerated |
| Predicted Species | Human |
| Human Gene ID | 11274 |
| Product Grade | https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png |
| ELISA Detection Limit | Antibody detection limit dilution 1:2000. |
| Western Blot | Western Blot: No signal obtained yet but low background observed in Human Liver extracts at up to 2µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? |