

GOAT ANTI-USP21 ANTIBODY

SKU: EB05998



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	ubiquitin-specific processing protease 23 ubiquitin carboxyl-terminal hydrolase 23 ubiquitin-specific protease 23 ubiquitin-specific protease 16 ubiquitin thiolesterase 23 deubiquitinating enzyme 23 NEDD8-specific protease MGC3394 USP23 USP16 ubiquitin specific protease 21 USP21
Accession ID	NP_036607; NP_057656
Immunogen	Peptide with sequence LFYQLMQEPPRCL, from the C Terminus of the protein sequence according to NP_036607; NP_057656.
Product Comments	This antibody is expected to recognise both human isoforms of this protein.
Peptide Sequence	LFYQLMQEPPRCL
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	27005
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:4000.
Western Blot	Western Blot: No signal obtained yet but low background observed in Human Heart and Muscle lysates at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA