

GOAT ANTI-USP21 ANTIBODY

SKU: EB08660



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	ubiquitin-specific protease 21 ubiquitin-specific protease 16 ubiquitin-specific processing protease 21 ubiquitin thiolesterase 21 ubiquitin specific protease 21 ubiquitin carboxyl-terminal hydrolase 21 deubiquitinating enzyme 21 NEDD8-specific protease USP23 USP16 MGC3394 ubiquitin specific peptidase 21 USP21
Accession ID	NP_036607.3
Blocking Peptide	EBP08660
Immunogen	Peptide with sequence GRTREPPVNIQP-C, from the N Terminus of the protein sequence according to NP_036607.3.
Product Comments	Reported variants represent identical protein: NP_001014443.1 NP_036607.3
Peptide Sequence	GRTREPPVNIQP-C
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Cow
Human Gene ID	27005
Mouse Gene ID	30941
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:64000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)