

GOAT ANTI-USP25 ANTIBODY

SKU: EB06002



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	ubiquitin specific protease USP25 USP21 ubiquitin specific protease 25 USP25
Accession ID	NP_037528.3
Immunogen	Peptide with sequence C-RIMLSLSRTPADGR, from the C Terminus of the protein sequence according to NP_037528.3.
Peptide Sequence	C-RIMLSLSRTPADGR
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat
Human Gene ID	29761
Mouse Gene ID	30940
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:4000.
Western Blot	Western Blot: No signal obtained yet but low background observed in Human Kidney and Testis extracts at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA