

# GOAT ANTI-USP26 ANTIBODY

**SKU:** EB06003



## SPECIFICATIONS

|                              |                                                                                                                                                                                                                                                                           |
|------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>             | 100 µg                                                                                                                                                                                                                                                                    |
| <b>Synonym / Alias Names</b> | ubiquitin-specific protease 26 ubiquitin-specific processing protease 26 ubiquitin thiolesterase 26 ubiquitin carboxyl-terminal hydrolase 26 deubiquitinating enzyme 26 OTTHUMP00000024046 MGC120068 MGC120067 MGC120066 RP1-97K10.1 ubiquitin specific protease 26 USP26 |
| <b>Accession ID</b>          | NP_114113.1                                                                                                                                                                                                                                                               |
| <b>Blocking Peptide</b>      | EBP06003                                                                                                                                                                                                                                                                  |
| <b>Immunogen</b>             | Peptide with sequence C-QLNSKEVEETLQKE, from the C Terminus of the protein sequence according to NP_114113.1.                                                                                                                                                             |
| <b>Peptide Sequence</b>      | C-QLNSKEVEETLQKE                                                                                                                                                                                                                                                          |
| <b>Shipping Instructions</b> | Refrigerated                                                                                                                                                                                                                                                              |
| <b>Predicted Species</b>     | Human                                                                                                                                                                                                                                                                     |
| <b>Human Gene ID</b>         | 83844                                                                                                                                                                                                                                                                     |
| <b>Product Grade</b>         | <a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png</a>                                   |
| <b>ELISA Detection Limit</b> | Antibody detection limit dilution 1:32000.                                                                                                                                                                                                                                |
| <b>Western Blot</b>          | Preliminary experiments gave no signal but low background in Human Testis lysates at up to 0.5µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?                                            |
| <b>Application Type</b>      | Pep-ELISA                                                                                                                                                                                                                                                                 |