

GOAT ANTI-USP28 (C-TERM) ANTIBODY

SKU: EB06004



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	ubiquitin specific peptidase 28 KIAA1515 ubiquitin specific protease 28 USP28
Accession ID	NP_065937.1
Blocking Peptide	EBP06004
Immunogen	Peptide with sequence C-MESIQGVSTVTVK, from the C Terminus of the protein sequence according to NP_065937.1.
Peptide Sequence	C-MESIQGVSTVTVK
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Pig, Cow
Human Gene ID	57646
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Western Blot	No signal obtained yet but low background observed in Human Testis and Skeletal Muscle lysates at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA