
Telephone: (650) 697-3600





GOAT ANTI-USP31 ANTIBODY

SKU: EB08665



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

Alias

ubiquitin specific proteinase 31|ubiquitin specific protease 31|KIAA1203|ubiquitin specific peptidase 31|USP31

Names

Accession

NP_065769.3

Blocking

EBP08665

Peptide

Peptide with sequence CKNTGDDEAERGHQ, from the internal region of the protein sequence according to

Immunogen NP 065769.3.

Peptide Sequence

CKNTGDDEAERGHQ

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Instructions

Refrigerated





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

Predicted

Human

Species

Reactive Human **Species**

Human Gene ID

57478

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

IHC Results

In paraffin embedded Human Cortex shows cellular body staining of neurons. Recommended concentration,

3-5μg/ml.

ELISA

Detection

Antibody detection limit dilution 1:8000.

Limit

Application

Pep-ELISA, IHC

Type

GALLERY IMAGES

