

GOAT ANTI-VASP (AA282-295) ANTIBODY

SKU: EB12618



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias Names VASP|vasodilator-stimulated phosphoprotein

Accession ID NP_003361.1

Blocking Peptide EBP12618

Immunogen Peptide with sequence C-EKTPKDESANQEPP, from the internal region of the protein sequence according to NP_003361.1.

Product Comments This antibody is expected to recognize.

Peptide Sequence C-EKTPKDESANQEPP

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species Human

Reactive Species Human

Human Gene ID 7408

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA Detection Limit Antibody detection limit dilution 1:1000.

Western Blot Approx 50kDa band observed in Human Platelets lysates (calculated MW of 39.8kDa according to NP_003361.1). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Reinhard et al, EMBO J. 1992 Jun;11(6):2063-70. PMID: 1318192). Recommended concentration: 0.1-0.3µg/ml.

Application Type Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

