

GOAT ANTI-VAV3 ANTIBODY

SKU: EB05250



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	vav 3 oncogene VAV3
Accession ID	NP_006104
Blocking Peptide	EBP05250
Immunogen	Peptide with sequence C-KQTEEKYTETLESI, from the C Terminus of the protein sequence according to NP_006104.
Peptide Sequence	C-KQTEEKYTETLESI
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, chicken
Reactive Species	Human
Human Gene ID	10451
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:4000.
Western Blot	Western Blot: Approx 97kDa band seen in Jurkat cell lysates. Recommended for use at 2-4µg/ml.