



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-VCAM1 ANTIBODY

SKU: EB12031



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / **Alias**

VCAM1|vascular cell adhesion molecule 1|CD106|INCAM-100|CD106 antigen|vascular cell adhesion protein 1

Names

Accession

NP_001069.1; NP_001186763.1

Blocking

Peptide

EBP12031

Immunogen

Peptide with sequence C-TQIDSPLSGKVRSE, from the internal region of the protein sequence according to

NP_001069.1; NP_001186763.1.

Product

This antibody is expected to recognize isoform a (NP 001069.1) and isoform c (NP 001186763.1).

Comments Peptide

Sequence

C-TQIDSPLSGKVRSE

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Human

Reactive

Human **Species**

Human

7412 Gene ID

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite medium.png

ELISA

Detection

Western

Antibody detection limit dilution 1:64000.

Limit

Approx 85kDa band observed in Human Uterus lysates (calculated MW of 81.3kDa according to NP_001069.1).

Recommended concentration: 1-3µg/ml. Blot

Application

Type

Pep-ELISA, WB

GALLERY IMAGES



Email: customerservice@vectorlabs.com

Telephone: <u>(650)</u> 697-3600

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa