

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GOAT ANTI-VDAC2 (INTERNAL) ANTIBODY**

**SKU:** EB07409



Telephone: (650) 697-3600



250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

10kDa



Telephone: (650) 697-3600

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 μg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

FLJ23841|RP11-375G3.1|voltage-dependent anion channel 2|VDAC2

**Alias** Names

**Accession** 

NP\_003366.2

**Blocking** 

EBP07409 **Peptide** 

Peptide with sequence C-DSAKSKLTRNN, from the internal region of the protein sequence according to **Immunogen** 

NP 003366.2.

**Peptide** 

C-DSAKSKLTRNN Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

Human, Mouse, Rat, Dog **Species** 

Reactive

Human **Species** 

Human

7417 Gene ID

Mouse

22334 **Gene ID** 

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png Grade

**ELISA** 

**Detection** Antibody detection limit dilution 1:64000.

Limit

Western Approx 32kDa band observed in Human Heart lysates (calculated MW of 31.6 kDa according to NP\_003366.2).

**Blot** Recommended concentration: 0.3-1µg/ml.

**Application** 

Pep-ELISA, WB Type

## **GALLERY IMAGES**



Telephone: (650) 697-3600



250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

10kDa