

GOAT ANTI-VDAC2 (INTERNAL) ANTIBODY

SKU: EB07409

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

10kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	FLJ23841 RP11-375G3.1 voltage-dependent anion channel 2 VDAC2
Accession ID	NP_003366.2
Blocking Peptide	EBP07409
Immunogen	Peptide with sequence C-DSAKSKLTRNN, from the internal region of the protein sequence according to NP_003366.2.
Peptide Sequence	C-DSAKSKLTRNN
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog
Reactive Species	Human
Human Gene ID	7417
Mouse Gene ID	22334
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Approx 32kDa band observed in Human Heart lysates (calculated MW of 31.6 kDa according to NP_003366.2). Recommended concentration: 0.3-1µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

