


Telephone: (650) 697-3600

## **GOAT ANTI-VE-CADHERIN ANTIBODY**



e erest



## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / **Alias** 

CDH5|cadherin 5, type 2 (vascular endothelium)|7B4|CD144|7B4 antigen|VE-cadherin|cadherin 5, type 2, VEcadherin (vascular epithelium)|cadherin-5|cd144 antigen|endothelial-specific cadherin|vascular endothelial

Names cadherin

Usage Summary <strong>Immunofluorescence:</strong> Strong expression of the protein seen in the cytoplasm of HeLa cells. Recommended concentration: 10µg/ml. <strong>Flow Cytometry:</strong> Flow cytometric analysis of

HeLa cells. Recommended concentration: 10ug/ml.

**Accession** 

NP\_001786.2

**Blocking Peptide** 

EBP12667









Peptide with sequence C-KQARAHGKSVPE, from the internal region of the protein sequence according to **Immunogen** NP 001786.2.

**Peptide** C-KQARAHGKSVPE **Sequence** 

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

**Shipping** Instructions

Refrigerated

**Predicted** 

Human, Pig, Cow **Species** 

Reactive

**Species** 

Human

Human

1003 Gene ID

**Product** Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png

**ELISA** 

**Detection** Antibody detection limit dilution 1:128000.

Limit

Application

Pep-ELISA, IF, FC **Type** 

## **GALLERY IMAGES**

