

# GOAT ANTI-VINCULIN ANTIBODY

SKU: EB09513



## SPECIFICATIONS

|                              |                                                                                                                                                                                                                                   |
|------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>             | 100 µg                                                                                                                                                                                                                            |
| <b>Synonym / Alias</b>       | metavinculin OTTHUMP00000043885 OTTHUMP00000019861 MVCL CMD1W vinculin VCL                                                                                                                                                        |
| <b>Names</b>                 |                                                                                                                                                                                                                                   |
| <b>Accession ID</b>          | NP_054706.1; NP_003364.1                                                                                                                                                                                                          |
| <b>Blocking Peptide</b>      | EBP09513                                                                                                                                                                                                                          |
| <b>Immunogen</b>             | Peptide with sequence C-RDYLDGSRGILSGT, from the internal regions (near N Terminus) of the protein sequence according to NP_054706.1; NP_003364.1.                                                                                |
| <b>Product Comments</b>      | This antibody is expected to recognize both reported isoforms (NP_054706.1; NP_003364.1).                                                                                                                                         |
| <b>Peptide Sequence</b>      | C-RDYLDGSRGILSGT                                                                                                                                                                                                                  |
| <b>Shipping Instructions</b> | Refrigerated                                                                                                                                                                                                                      |
| <b>Predicted Species</b>     | Human, Mouse, Rat, Dog, Cow, Pig                                                                                                                                                                                                  |
| <b>Reactive Species</b>      | Mouse                                                                                                                                                                                                                             |
| <b>Human Gene ID</b>         | 7414                                                                                                                                                                                                                              |
| <b>Product Grade</b>         | <a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png</a> |
| <b>ELISA Detection Limit</b> | Antibody detection limit dilution 1:16000.                                                                                                                                                                                        |
| <b>Application Type</b>      | Pep-ELISA, WB                                                                                                                                                                                                                     |