

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-VN1R1 ANTIBODY

SKU: EB10484

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

VN1R1|ZVNR1|ZVNH1|vomeronasal olfactory receptor, (chromosome 19) subtype I, member 1|vomeronasal 1 **Alias**

receptor 1|VNR19I1|V3R-related|V1R-like 1|V1RL1 **Names**

Accession

NP 065684.1 ID

Blocking

EBP10484 **Peptide**

Peptide with sequence C-SSDLNENQHP, from the internal region of the protein sequence according to Immunogen

NP 065684.1.

Peptide

C-SSDLNENQHP Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Human

Human Gene ID

57191

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

Grade **ELISA**

Detection

Antibody detection limit dilution 1:16000.

Limit

Application Pep-ELISA **Type**