

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-VN1R1 ANTIBODY

SKU: EB10483

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

VN1R1|ZVNR1|ZVNH1|vomeronasal olfactory receptor, (chromosome 19) subtype I, member 1|vomeronasal 1 Alias

receptor 1|VNR19I1|V3R-related|V1R-like 1|V1RL1

Names

Accession NP 065684.1 ID

Blocking

EBP10483 **Peptide**

Immunogen Peptide with sequence QHNHSNRLSC, from the internal region of the protein sequence according to NP 065684.1.

Peptide

QHNHSNRLSC Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping Instructions

Refrigerated

Predicted

Human **Species**

Human

57191 Gene ID

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

ELISA

Detection

Antibody detection limit dilution 1:1000.

Limit

Application Pep-ELISA **Type**