

GOAT ANTI-VN1R1 ANTIBODY

SKU: EB10483



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	VN1R1 ZVNR1 ZVNH1 vomeronasal olfactory receptor, (chromosome 19) subtype I, member 1 vomeronasal 1 receptor 1 VNR19I1 V3R-related V1R-like 1 V1RL1
Accession ID	NP_065684.1
Blocking Peptide	EBP10483
Immunogen	Peptide with sequence QHNHSNRLSC, from the internal region of the protein sequence according to NP_065684.1.
Peptide Sequence	QHNHSNRLSC
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	57191
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:1000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)

