

GOAT ANTI-VTCN1 (AA153-165) ANTIBODY

SKU: EB11285

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	VTCN1 V-set domain-containing T-cell activation inhibitor 1 V-set domain containing T cell activation inhibitor 1 VCTN1 T-cell costimulatory molecule B7x T cell costimulatory molecule B7x PRO1291 immune costimulatory protein B7-H4 FLJ22418 B7X B7S1 B7-H4 B7H4 B7h.5
Accession ID	NP_078902.2
Blocking Peptide	EBP11285
Immunogen	Peptide with sequence C-PEVNVDYNASSET, from the internal region of the protein sequence according to NP_078902.2.
Peptide Sequence	C-PEVNVDYNASSET
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow
Reactive Species	Human
Human Gene ID	79679
Mouse Gene ID	242122
Rat Gene ID	295322
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Approx 28kDa band observed in Human Spleen lysates (calculated MW of 30.9kDa according to NP_078902.2). Recommended concentration: 0.1-0.3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

