


GOAT ANTI-WBP1 ANTIBODY

SKU: EB10322



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	WBP1 WW domain binding protein 1 WBP-1 MGC15305
Names	
Accession ID	NP_036609.1
Blocking Peptide	EBP10322
Immunogen	Peptide with sequence C-PAHFEGTNVE, from the internal region of the protein sequence according to NP_036609.1.
Peptide Sequence	C-PAHFEGTNVE
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	23559

Mouse Gene ID	22377
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Brain (Cerebellum) and Pancreas. Recommended concentration: 3.75µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:16000.
Application Type	Pep-ELISA, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

