


Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-WBP1 ANTIBODY

SKU: EB10322



SPECIFICATIONS

Unit Size 100 μg

Synonym /

Alias WBP1|WW domain binding protein 1|WBP-1|MGC15305

Names

Accession

NP_036609.1 ID

Blocking

EBP10322 **Peptide**

Peptide with sequence C-PAHFEGTNVE, from the internal region of the protein sequence according to

Peptide

C-PAHFEGTNVE Sequence

Shipping

Refrigerated Refrigerated

Predicted Species

Human

Reactive

Human

Species Human

23559 **Gene ID**





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

Mouse **Gene ID**

22377

Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. a mazon aws.com/assets/products/image/elite_medium.png$ Grade

IHC Results Paraffin embedded Human Brain (Cerebellum) and Pancreas. Recommended concentration: 3.75µg/ml.

ELISA

Detection

Antibody detection limit dilution 1:16000.

Limit

Application Pep-ELISA, IHC

Type

GALLERY IMAGES

