

GOAT ANTI-WBP1 ANTIBODY

SKU: EB10321



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	WBP1 WW domain binding protein 1 WBP-1 MGC15305
Names	
Accession ID	NP_036609.1
Blocking Peptide	EBP10321
Immunogen	Peptide with sequence C-PFPTGSLLDLR, from the internal region of the protein sequence according to NP_036609.1.
Peptide Sequence	C-PFPTGSLLDLR
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow
Human Gene ID	23559
Mouse Gene ID	22377
Rat Gene ID	297381
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Preliminary experiments gave an approx 40kDa band in Human Pancreas lysates after 1µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 29.1kDa according to NP_036609.1 The 40kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)