

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-WBP2NL / PAWP ANTIBODY

SKU: EB07830

SPECIFICATIONS

Unit Size

 $100 \mu g$

Synonym / Alias

postacrosomal sheath WW domain-binding

Names

protein|OTTHUMP00000028745|PAWP|MGC26816|FLJ26145|CTA-250D10.11|WBP2 N-terminal like|WBP2NL

Accession

ID

NP_689826.1

Blocking

Peptide

EBP07830

Immunogen

Peptide with sequence NQSHTENRRGALIPN-C, from the N Terminus of the protein sequence according to

NP 689826.1.

Peptide

Sequence

NQSHTENRRGALIPN-C

Shipping

Refrigerated Instructions

Human

Predicted

Species

Human

Gene ID

164684

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

ELISA

Detection

Antibody detection limit dilution 1:16000.

Limit

Western Blot

Preliminary experiments in Human Testis lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been

reported with other antibodies/lysates?

Application

Type

Pep-ELISA