

GOAT ANTI-WBP2NL / PAWP ANTIBODY

SKU: EB07830



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	postacrosomal sheath WW domain-binding protein OTTHUMP00000028745 PAWP MGC26816 FLJ26145 CTA-250D10.11 WBP2 N-terminal like WBP2NL
Accession ID	NP_689826.1
Blocking Peptide	EBP07830
Immunogen	Peptide with sequence NQSHTENRRGALIPN-C, from the N Terminus of the protein sequence according to NP_689826.1.
Peptide Sequence	NQSHTENRRGALIPN-C
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	164684
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Preliminary experiments in Human Testis lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)

