Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-WDR9 / BRWD1 (AA1441-1455) ANTIBODY**

**SKU:** EB11490

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% boyine serum albumin.

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Alias

Synonym / WDR9|WD repeat-containing protein 9|WD repeat protein WDR9-form2|transcriptional unit

Names and WD repeat-containing protein 1|bromodomain and WD repeat domain containing 1  $\,$ 

Accession

NP\_061836.2; NP\_387505.1; NP\_001007247.1

**Blocking** 

ID

EBP11490 Peptide

Peptide with sequence RQNCKGDSQPNKSIR, from the internal region of the protein sequence according to NP\_061836.2; NP\_387505.1; Immunogen

NP\_001007247.1.

**Product** Comments

This antibody is expected to recognize reported isoforms A and B (NP 061836.2; NP 387505.1).

Peptide

RQNCKGDSQPNKSIR Sequence

Purification

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide. Method

**Shipping** 

Refrigerated Instructions

Predicted

Human Species

Human

54014

Gene ID

Product Grade

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png\\$ 

Detection

Limit

Type

Antibody detection limit dilution 1:2000

**Application** Pep-ELISA