

GOAT ANTI-WDR9 / BRWD1 (AA1441-1455) ANTIBODY

SKU: EB11490



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / WDR9|WD repeat-containing protein 9|WD repeat protein WDR9-form2|transcriptional unit

Alias N143|OTTHUMP00000207967|OTTHUMP00000068928|OTTHUMP00000068797|OTTHUMP00000068796|N143|FLJ43918|C21orf107|BRWD1|bromodomain and WD repeat-containing protein 1|bromodomain and WD repeat domain containing 1

Accession ID NP_061836.2; NP_387505.1; NP_001007247.1

Blocking Peptide EBP11490

Immunogen Peptide with sequence RQNCKGDSQPNSIR, from the internal region of the protein sequence according to NP_061836.2; NP_387505.1; NP_001007247.1.

Product Comments This antibody is expected to recognize reported isoforms A and B (NP_061836.2; NP_387505.1).

Peptide Sequence RQNCKGDSQPNSIR

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species Human

Human Gene ID 54014

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

ELISA Detection Limit Antibody detection limit dilution 1:2000.

Application Type Pep-ELISA

DOCUMENTS

- [Data Sheet](#)