

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600


GOAT ANTI-WNT15 / WNT9B ANTIBODY

SKU: EB09452



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600







Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

wingless-type MMTV integration site family, member 15|WNT15|WNT14B|wingless-type MMTV integration site **Alias**

family, member 9B|WNT9B Names

Accession

NP_003387.1 ID

Blocking

EBP09452 **Peptide**

Peptide with sequence C-KRGNKDLRARADA, from the internal region of the protein sequence according to **Immunogen**

NP 003387.1.

Peptide Sequence

C-KRGNKDLRARADA

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Mouse, Rat **Species**

Reactive

Human, Mouse **Species**

Human

7484 Gene ID

Mouse

22412 **Gene ID**

Rat Gene ID 303586

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Grade ELISA

Detection Antibody detection limit dilution 1:8000.

Limit

Western Approx 38kDa band observed in Human Brain Cerebellum lysates and in Mouse Kidney and Brain lysates

(calculated MW of 39kDa according to NP_003387.1). Recommended concentration: 0.3-1µg/ml. Blot

Application

Pep-ELISA, WB **Type**

> Goat Anti-WNT15 / WNT9B Antibody https://everestbiotech.com/products/goat-anti-wnt15-wnt9b-antibody/



Telephone: (650) 697-3600



GALLERY IMAGES

