

GOAT ANTI-WNT2 ANTIBODY

SKU: EB05372



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	secreted growth factor INT1L1 IRP wingless-type MMTV integration site family member 2 WNT2
Accession ID	NP_003382
Blocking Peptide	EBP05372
Immunogen	Peptide with sequence SWWYMRATGGSSRV-C, from the N Terminus of the protein sequence according to NP_003382.
Peptide Sequence	SWWYMRATGGSSRV-C
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat
Human Gene ID	7472
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
Western Blot	Western blot: No signal obtained yet but low background observed in Human Brain and NIH-3T3 cell lysates at upto 2µg/ml.