Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-XAGE1 (ISOFORMS 1A, 1B AND 1C) ANTIBODY

SKU: EB06435

SPECIFICATIONS

Unit Size $100 \mu g$

Synonym /

Alias XAGE-1 protein|G antigen, family D, 2|XAGE-1|XAGE1|GAGED2

Names

Accession

NP_065144; CAC38108

Blocking Peptide

EBP06435

Peptide with sequence CPEAGEEQPQV, from the C Terminus of the protein sequence according to NP_065144; **Immunogen**

CAC38108.

Refrigerated

Product

NCBI Reference Sequence: NP_065144.2 removed as redundant sequence Comments

Peptide

CPEAGEEQPQV

Sequence

Shipping

Instructions

Species

Predicted Human

Human

9503

Gene ID Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

Grade **ELISA**

Detection Antibody detection limit dilution 1:2000.

Blot

Limit

Preliminary experiments in Human Testis lysates and in lysates of different Human epithelial cell lines gave no specific signal but low background (at antibody concentration up to 1µg/ml).

Application

Type

Western

Pep-ELISA