

# GOAT ANTI-XAGE1 (ISOFORMS 1A, 1B AND 1C) ANTIBODY

SKU: EB06435



## SPECIFICATIONS

<b>Unit Size</b>	100 µg
<b>Synonym /</b>	
<b>Alias</b>	XAGE-1 protein G antigen, family D, 2 XAGE-1 XAGE1 GAGED2
<b>Names</b>	
<b>Accession ID</b>	NP_065144; CAC38108
<b>Blocking Peptide</b>	EBP06435
<b>Immunogen</b>	Peptide with sequence CPEAGEEQPV, from the C Terminus of the protein sequence according to NP_065144; CAC38108.
<b>Product Comments</b>	NCBI Reference Sequence: NP_065144.2 removed as redundant sequence
<b>Peptide Sequence</b>	CPEAGEEQPV
<b>Shipping Instructions</b>	Refrigerated
<b>Predicted Species</b>	Human
<b>Human Gene ID</b>	9503
<b>Product Grade</b>	<a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png</a>
<b>ELISA</b>	
<b>Detection Limit</b>	Antibody detection limit dilution 1:2000.
<b>Western Blot</b>	Preliminary experiments in Human Testis lysates and in lysates of different Human epithelial cell lines gave no specific signal but low background (at antibody concentration up to 1µg/ml).
<b>Application Type</b>	Pep-ELISA

## DOCUMENTS

- [Data Sheet](#)