

GOAT ANTI-YAP / YAP1 ANTIBODY

SKU: EB09727



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	YAP1 Yes-associated protein 1, 65kDa YAP YAP2 YAP65 YKI yes-associated protein 2
Names	
Accession ID	NP_001123617.1; NP_006097.1
Immunogen	Peptide with sequence C-NQSTLPSQQNRFPD, from the internal region of the protein sequence according to NP_001123617.1; NP_006097.1.
Product Comments	This antibody is expected to recognize both reported isoforms (NP_001123617.1; NP_006097.1).
Peptide Sequence	C-NQSTLPSQQNRFPD
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow
Reactive Species	Rat
Human Gene ID	10413
Mouse Gene ID	22601
Rat Gene ID	363014
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Approx 55kDa band observed in Rat Muscle lysates (calculated MW of 54.5kDa according to NP_001123617.1). Recommended concentration: 0.1-0.5µg/ml.
Application Type	ELISA, WB

SELECTED REFERENCES

[{"pmid": 18953429, "intro": "", "title": "The neurofibromatosis 2 tumor suppressor gene product, merlin, regulates human meningioma cell growth by signaling through YAP.", "author": "Striedinger K, VandenBerg SR, Baia GS, McDermott MW, Gutmann DH, Lal A.", "journal": "Neoplasia. 2008 Nov;10(11):1204-12."}]

GALLERY IMAGES

