


GOAT ANTI-ZBTB7B ANTIBODY

SKU: EB08598



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names zinc finger protein 67 homolog|zinc finger and BTB domain containing 15|kruppel-related zinc finger protein hcKrox|hcKrox|c-Krox|ZFP67|ZBTB15|DKFZp686G01254|zinc finger and BTB domain containing 7B|ZBTB7B

Accession ID NP_056956.2

Blocking Peptide EBP08598

Immunogen Peptide with sequence C-DSQDKLVRKRRSQ, from the internal region of the protein sequence according to NP_056956.2.

Peptide Sequence C-DSQDKLVRKRRSQ

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species	Human, Mouse, Rat, Dog, Cow
Reactive Species	Human
Human Gene ID	51043
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	In paraffin embedded Human Testis shows strong nuclear staining in primary spermatocytes. Recommended concentration, 3-5µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:64000.
Application Type	Pep-ELISA, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

