


Telephone: (650) 697-3600

GOAT ANTI-ZCCHC11 / PAPD3 ANTIBODY

SKU: EB05042



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage Instructions

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias PAP associated domain containing 3|PAPD3|zinc finger, CCHC domain containing 11|KIAA0191|ZCCHC11

Names

Usage Immunofluorescence: Strong expression of the protein seen in the cytoplasm of A431 cells.

Summary Recommended concentration: 10µg/ml.

Accession

NP_001009881.1; NP_056084.1; NP_001009882.1

Blocking

EBP05042 **Peptide**

Peptide with sequence RERCPHPPRGNVSE, from the C Terminus of the protein sequence according to **Immunogen**

NP 001009881.1; NP 056084.1; NP 001009882.1.









Product This antibody is expected to recognize all 3 reported isoforms (as represented by NP_001009881.1;

Comments NP_056084.1; NP_001009882.1)

Peptide Sequence

RERCPHPPRGNVSE

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping Instructions

Refrigerated

Predicted

Species Human, Mouse, Rat, Dog, Cow

Reactive

Species

Human

Human

Gene ID

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

IHC Results

In paraffin embedded Human Lymph Node shows nuclear staining of select lymphocytes and staining of

macrophages. Recommended concentration: 2-6µg/ml.

ELISA

Detection

Antibody detection limit dilution 1:16000.

Limit

Application

Туре

Pep-ELISA, IF, IHC

GALLERY IMAGES

