

GOAT ANTI-ZIC4 ANTIBODY

SKU: EB05746



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	zinc finger protein of the cerebellum 4 zinc finger protein ZIC 4 zinc family member 4 protein HZIC4 ZIC4
Accession ID	AF332509_1
Blocking Peptide	EBP05746
Immunogen	Peptide with sequence C-EAGPRVWISSPN, from the C Terminus of the protein sequence according to AF332509_1.
Product Comments	The immunizing peptide no longer represents the C-terminal sequence of the original reported isoform NP_115529.
Peptide Sequence	C-EAGPRVWISSPN
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	84107
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Western Blot	Preliminary experiments in Human Brain (Cerebellum) lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA