


Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-ZIC5 ANTIBODY

SKU: EB05745



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Instructions

Synonym / **Alias**

Zic family member 5 (odd-paired homolog, Drosophila)|zinc family member 5 protein|zinc finger protein of the

cerebellum 5|ZIC5 **Names**

Accession

NP_149123.2

Blocking

ID

EBP05745

Peptide

Immunogen

Peptide with sequence C-EIYGNPEVVRTIH, from the C Terminus of the protein sequence according to NP_149123.2.

Peptide

C-EIYGNPEVVRTIH

Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

Predicted

Human

Species Reactive

Human **Species**

Human

Gene ID

85416

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

IHC Results Paraffin embedded Human Small Intestine, Skin and Prostate. Recommended concentration: 3.75µg/ml.

ELISA

Detection

Antibody detection limit dilution 1:128000.

Limit

Type

Application

Pep-ELISA, IHC

GALLERY IMAGES

