

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-ZNF228 / ZFP112 ANTIBODY

SKU: EB09120

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias zinc finger protein 228|ZNF228|zinc finger protein 112 homolog (mouse)|ZFP112

Names

Accession NP 001076804.1; NP 037512.3 ID

Blocking

EBP09120 **Peptide**

Peptide with sequence C-EREDDIENSHLKSYQ, from the internal region of the protein sequence according to Immunogen

NP 001076804.1; NP 037512.3.

Product Comments

This antibody is expected to recognize both reported isoforms (NP_001076804.1; NP_037512.3)

Peptide

C-EREDDIENSHLKSYQ Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human Species

Human

7771

Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:4000.

Limit

Application Pep-ELISA **Type**