


GOAT ANTI-ZNF230 / RNF141 ANTIBODY

SKU: EB05954



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias RNF141|ZNF230|ZFP26|ring finger protein 141|C3HC4-like zinc finger protein|MGC8715

Names

Accession ID NP_057506.2

Blocking Peptide EBP05954

Immunogen Peptide with sequence C-LNMADEAGQPHRP, from the C Terminus of the protein sequence according to NP_057506.2.

Peptide Sequence C-LNMADEAGQPHRP

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species Human, Mouse, Rat, Dog, Cow

Reactive Species Human

Human Gene ID 50862

Mouse Gene ID 67150

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA Detection Limit Antibody detection limit dilution 1:16000.

Western Blot Approx 25-30kDa band observed in Human Testis lysate (predicted MW of 26kDa according to NP_057506). Recommended for use at 1µg/ml.

Application Type Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

