



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-ZNF703 ANTIBODY

SKU: EB06895







250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias

ZNF703|zinc finger protein 703|FLJ14299|ZNF503L|hypothetical protein LOC80139

Names

ID

Accession

NP_079345.1

Blocking

EBP06895 **Peptide**

Peptide with sequence C-ELDKKDQEPKPSPE, from the internal region of the protein sequence according to **Immunogen**

NP 079345.1.

Peptide

C-ELDKKDQEPKPSPE Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Dog, Pig **Species**

Reactive

Human **Species**

Human

80139 **Gene ID**

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:32000.

Limit

Western Approx 60kDa band observed in Human Breast cancer lysates (calculated MW of 58.2kDa according to

Blot NP_079345.1). Recommended concentration: 1-3µg/ml.

Application

Pep-ELISA, WB Type

GALLERY IMAGES



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa