


# GOAT ANTI-ZNF703 ANTIBODY

**SKU:** EB06895



## SPECIFICATIONS

**Formulation** Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 µg

**Storage Instructions** Aliquot and store at -20°C. Minimize freezing and thawing.

**Synonym /**

**Alias** ZNF703|zinc finger protein 703|FLJ14299|ZNF503L|hypothetical protein LOC80139

**Names**

**Accession ID** NP\_079345.1

**Blocking Peptide** EBP06895

**Immunogen** Peptide with sequence C-ELDKKDQEPKPSPE, from the internal region of the protein sequence according to NP\_079345.1.

**Peptide Sequence** C-ELDKKDQEPKPSPE

**Purification Method** Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

**Shipping Instructions** Refrigerated

**Predicted Species** Human, Dog, Pig

**Reactive Species** Human

**Human Gene ID** 80139

**Product Grade** [https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\\_medium.png](https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png)

**ELISA**

**Detection Limit** Antibody detection limit dilution 1:32000.

**Western Blot** Approx 60kDa band observed in Human Breast cancer lysates (calculated MW of 58.2kDa according to NP\_079345.1). Recommended concentration: 1-3µg/ml.

**Application Type** Pep-ELISA, WB

## DOCUMENTS

- [Data Sheet](#)

## GALLERY IMAGES

