


MOUSE ANTI-GFP MONOCLONAL ANTIBODY 168AT1211

SKU: EB00102



SPECIFICATIONS

Formulation Purified monoclonal antibody is supplied in PBS with 0.09% (W/V) sodium azide.

Unit Size 100 µg

Storage Instructions Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C.

Synonym /

Alias Green Fluorescent Protein

Names

Usage Summary The suggested dilutions: ELISA: 1:3000 Western blotting: 1:2000 Immunofluorescence: 1:50

Immunogen This antibody is generated from His-tagged GFP protein.

Purification Method This antibody is purified through a protein G column, eluted with high and low pH buffers and neutralized immediately, followed by dialysis against PBS.

Shipping Instructions Refrigerated

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Application
Type ELISA, WB, IF

GALLERY IMAGES

