


Telephone: (650) 697-3600

MOUSE ANTI-HIS TAG MONOCLONAL ANTIBODY 6AT18

SKU: EB00105





Telephone: (650) 697-3600

SPECIFICATIONS

Unit Size 100 μg

Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C.

Usage

The suggested dilutions: ELISA: 1:5000 Western blotting: 1:1000 **Summary**

This antibody was generated from a purified recombinant HIS-tagged fusion protein and poly-HIS **Immunogen**

(hexahistidine) peptide.

Purification

Purified monoclonal antibody is supplied in PBS with 0.09% (W/V) sodium azide.

Shipping

Refrigerated Instructions

Product Grade

Method

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Application

ELISA, WB **Type**

GALLERY IMAGES





Telephone: (650) 697-3600



