


MOUSE ANTI-HIS TAG MONOCLONAL ANTIBODY 6AT18

SKU: EB00105



SPECIFICATIONS

Unit Size	100 µg
Storage Instructions	Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C.
Usage Summary	The suggested dilutions: ELISA: 1:5000 Western blotting: 1:1000
Immunogen	This antibody was generated from a purified recombinant HIS-tagged fusion protein and poly-HIS (hexahistidine) peptide.
Purification Method	Purified monoclonal antibody is supplied in PBS with 0.09% (W/V) sodium azide.
Shipping Instructions	Refrigerated
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Application Type	ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

