

POSITIVE CONTROL LYSATE FOR ANTI- AKT3 ANTIBODY

SKU: EBL020-AKT3-P2



SPECIFICATIONS

Formulation	5ug/ul protein in RIPA buffer
Unit Size	50 µl
Storage Instructions	Store lysate at -20°C upon arrival. Avoid repeated freeze /thawed cycles.
Usage Summary	Western Blot: 1µg/ml staining of Primary antibody in Nuclear lysates EBL020-AKT3-P2, approx.55kDa band observed. (Calculated MW of 55.8kDa according to NP_005456.1). Primary incubation 1 hour at room temperature. Secondary antibody is EB2ND-002-A at 1:40K.
Shelf Life	3 months
Shipping Instructions	Refrigerated
Application Type	WB

GALLERY IMAGES

