

# **POSITIVE CONTROL LYSATE FOR ANTI- AKT3 ANTIBODY**

**SKU:** EBL020-AKT3-P1



---

## SPECIFICATIONS

<b>Formulation</b>	5ug/ul protein in RIPA buffer
<b>Unit Size</b>	50 µl
<b>Storage Instructions</b>	Store lysate at -20°C upon arrival. Avoid repeated freeze /thawed cycles.
<b>Usage Summary</b>	Western Blot: 1µg/ml staining of Primary antibody in Nuclear lysates EBL020-AKT3-P1, approx.50-55kDa band observed. (Calculated MW of 54kDa according to Human NP_859029.1). Primary incubation 1 hour at room temperature. Secondary antibody is EB2ND-002-A at 1:40K.
<b>Shelf Life</b>	3 months
<b>Shipping Instructions</b>	Refrigerated
<b>Application Type</b>	WB

## GALLERY IMAGES

